

# PORTABLE OUTLET

**Model: PO-160UPS**  
Uninterruptible Power Supply  
and CPAP Battery


## 159W Rechargeable UPS and CPAP Battery

- Works with the most popular CPAPs and BiPAPs
- Provides 1.5 to 2 nights of therapy on average, without heated tube or heated humidifier
- Does not require additional connectors or adapters
- Pure Sine Wave Technology adds battery life
- Optional Solar Panel recharges in a few hours
- Color Display & Power Percentage Indicator
- Meets TSA & FAA guidelines for air travel
- Can power laptops and other electronics
- Fast Charging USB - Output 5V/3.0 A
- AC Discharge Output: 110V/60Hz
- 1 year warranty




***When power goes out, your CPAP stays ON!***

Plug your CPAP into the Portable Outlet.  
The battery stays fully charged.



**When there is an outage, the battery switches ON.**  
Your therapy continues without any interruption.



**Also a great stand-alone battery for camping, travel or emergencies!**

Optional Solar Panel recharges in a few hours.



Includes travel bag and power cord.



Can power a laptop & phone.



Visit [www.PortableOutlet.com](http://www.PortableOutlet.com) or call 470-297-5778 to learn more.